# AEG


## PWI 4914 FA

### Bedienungsanleitung/Garantie

Gebbruksaanwijzing/Garantie • Mode d'emploi/Garantie  
Istruzioni per l'uso/Garanzia • Instrucciones de servicio/Garantía  
Manual de instruções/Garantia • Instruction Manual/Guarantee  
Bruksanvisning/Garanti • Instrukcja obsługi/Gwarancja  
Návod k použití/Záruka • A használati utasítás/garancia  
Руководство по эксплуатации/Гарантия

### Infrarot-Körperanalysewaage

Infrarode lichaamsanalyseweegschaal •  
Balance impédancemètre à infrarouge •  
Bilancia per l'analisi del corpo a raggi infrarossi •  
Báscula impedancímetro con infrarrojos •  
Balança infravermelha para análise de gordura do corpo •  
Infrared Body Analysis Scale •  
Infrarød kroppsanalysevekt • Waga na podczerwień z funkcją analizy •  
Infračervená osobní váha s analyzátozem tělesného tuku •  
Infravörös testelemző mérleg •  
Инфракрасные весы анализатор организма

CE

# Inhalt

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## General Safety Instructions

- Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing.
- The appliance is designed exclusively for private use and for the envisaged purpose. This appliance is not fit for commercial use. Do not use it outdoors (except if it is designed to be used outdoors). Keep it away from sources of heat, direct sunlight, humidity (never dip it into any liquid) and sharp edges. Do not use the appliance with wet hands.
- Check the appliance and the cable for damage on a regular basis. Do not use the appliance if it is damaged.
- Do not try to repair the appliance on your own. Always contact an authorized technician.
- Pay careful attention to the following "Special Safety Instructions".

## Children and Frail Individuals

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach. **Caution!** Do not allow small children to play with the foil as there is a **danger of suffocation!**
- In order to protect children or frail individuals from the hazards of electrical equipment, please ensure that this device is only used under supervision. This device is not a toy. Do **not** allow small children to play with it.

## Special safety notes

- The body fat scales are not suitable for the following people:
  - Children under the age of 6 years
  - Individuals with symptoms of fever, oedemas or osteoporosis
  - People undergoing dialysis treatment
  - Those taking cardiovascular medication
  - Pregnant women

- People with a pulse of less than 60 beats a minute
- The device can bear a maximum body weight of 150 kg (330,7 lb.)! A larger weight can damage the device!

## Overview of the Components

- 1 On switch **MODE**
- 2 **P/▲ / ATH/▼** Selection arrow keys
- 3 **SET** Storage button for personal values
- 4 **UNIT** button (setting)
- 5 LCD Display
- 6 Contact surfaces (display)
- 7 Infrared lamp
- 8 Indicator lamp

### Rear side

- 9 Battery Compartment

### Display

- 10 Storage location/date
- 11 Age (**AGE**)
- 12 **HT** body height cm / in(ch)
- 13 Percentage – body fat / body water / muscle mass
- 14 **WT** body weight kg, lb. St:
- 15 Calorie requirement **KCAL**
- 16 Body fat display **BF**
- 17 Body water display **BW**
- 18 Indication of muscle mass **BM**
- 19 Indication of **BON** mass
- 20 Male
- 21 Female
- 22 Child
- 23 Symbol 'Athlete'

## Notes on Handling Batteries

- Keep batteries away from children. They are **not** a toy!
- Never throw batteries into a fire or immerse in water.
- When inserting a battery, ensure the correct polarity.
- Do not open the batteries by force.
- Avoid contact with metallic objects. (Rings, nails, screws et cetera.) There is a **danger of short-circuiting!**
- As a result of a short-circuit batteries may heat up considerably or even catch fire. This may result in burns.
- For your own safety the battery poles

should be covered with adhesive strips during transport.

- If a battery leaks do not rub the liquid into your eyes or mucous membranes. If hands come into contact with the liquid they should be washed thoroughly and eyes rinsed out with clear water; if symptoms persist please consult a doctor.
- Used batteries should not be disposed of in the normal domestic waste. Please return old batteries to an appropriate connection point.

Germany:

- As a consumer you are required by law to return used (rechargeable) batteries.
- You can return your old (rechargeable) batteries free of charge where they were bought, as well as at the local public collection points in your town.

The following symbols are to be found on (rechargeable) batteries containing harmful substances:

Pb = contains lead  
Cd = contains cadmium  
Hg = contains mercury  
Li = contains lithium

## Installing the batteries (batteries not included)

1. Open the battery compartment on the rear side of the scale and the LCD.
2. Put in 3 batteries of type AM4 "AAA" 1.5V. Please ensure that the polarity is correct (see the floor of the battery compartment)!
3. Close the battery compartment again.
4. If the device is not used for a prolonged period of time, please take the batteries out in order to prevent leakage of battery acid.
5. Different battery types or new and used batteries must not be used together.
6. If the message "Lo-1" (LCD) or "Lo-2" (scale) for low battery should appear on the display replace the batteries.

### Attention!

Batteries do not belong in your usual household rubbish. Please dispose of old batteries at proper collecting points or at your dealer.

### Warning!

Do not expose batteries to high temperatures or direct sunlight. Never throw batteries into a fire.

**Danger of explosion!**

## Wall mounting

**Attention:** The LCD display is suitable for wall mounting. Use the included mounting material to attach the LCD display to a wall. Mount the display in a place where you can conveniently and safely operate it. First, mount the mounting plate; use the plate as a template to mark the mounting points. Align the mounting plate using a spirit level and mark the positions of the three holes on the wall. Drill three holes in the marked spots. Make sure not to damage any utility lines running through the wall! Hook the display (2 hooks) onto the mounting plate.

## Information for body fat measurement

People who use a heart pacemaker or another type of medical implant with electronic components are advised not to use these scales, as otherwise the performance of the implant may be impaired or endangered.

The device functions correctly if your feet are bare and dry and if you are in standing position.

The basic principle of these scales is measurement of the electrical impedance of the body. The scales send a weak, imperceptible electrical signal through your body in order to determine the current values. Under normal circumstances the signal is completely harmless! Nevertheless you should read our special safety instructions! The measurement is known as "biometric impedance analysis"

## GB

(BIA), a factor which is associated with the body fat and body water ratio and which is subject to other biological aspects (age, sex, height).

There are heavy people with a well-balanced relation between fat and muscular mass and there are apparently slender people with substantially higher body fat values. Fat is a factor that contributes to burdening the body and, in particular, the circulatory system. The checks of both values (fat and weight) are hence a basic prerequisite for health and fitness.

### Body Fat

Body fat is an important part of the body. Fat plays an important role in your body ranging from the protection of joints with its reserve of vitamins to the regulation of the body temperature. The goal is not to reduce drastically the quantity of fat in your body, but to reach and maintain a balanced relation between the muscular mass and the fat mass.

### The significance of water for our health

The human organism consists of approximately 55-60% water, with this quantity varying according to age and sex. Water has a number of tasks to fulfil in the body:

It is a building block of our cells. This means that all cells in the body, whether these are skin, gland, muscle, brain or other cells, can only function if they receive enough water. It is a solvent for important substances in our organism.

It is a means of transport for nutrients, the body's own substances and metabolic products. Please see the table of guideline values for the reference value for your fitness profile. It shows that there is a relation between the body fat values and your age and sex. Always bear in mind that, as in the case of the water balance in your body, also the share of fat is subject to great variations. Therefore, weigh yourself always at the same hour of the day to obtain a comparable value. The results of the analysis of the body fat and body water ratio are only reference values. They are not recommended for use for medical purposes. You should always consult your doctor if you have any questions about medical treatment or a diet to achieve your normal weight.

## Technical Features

Memory presets:

Memory	Sex	Height	Age(AGE)
01	Male	170 cm (67.0 in (inch))	30

Memory locations:	16
Height setting:	from 80 to 240 cm (31.5 – 94,5 in)
Age:	from 6 to 100 years
Minimum weight:	2.0 kg (4,4lb. / 0,314 St.)
Maximum weight:	150 kg (330,7lb. / 23,62 St.)
Weight measurement acc.:	100 g steps
Body fat gradation:	0.1%
Body fat measurement range:	1,0%–60,0%
Body water measurement range:	20,0%–75,0%
Muscle mass measuring range:	10,0%–50,0%
Calories measuring range:	0-9999 Kcal.
Calories gradation:	1 Kcal.
Time gradation:	24 hours
Permissible storage temperature:	–10°C to +60°C (14.0°F to 140.0°F)
Calendar:	2000–2099
Transmission distance:	max. 2m
Weight measuring accuracy:	+/-500g (+/-11lb)

Guidelines – body fat percentage

Age (years)	Women				Men			
	Under-weight	Normal	Slightly over-weight	Over-weight	Under-weight	Normal	Slightly over-weight	Over-weight
6-11	11-13%	14-23%	24-27%	>27%	9-11%	11-20%	21-24%	>24%
12-15	12-15%	16-24%	24-29%	>29%	7-9%	10-19%	20-23%	>23%
16-39	17-20%	21-28%	29-33%	>33%	10-12%	13-20%	21-25%	>25%
40-60	20-24%	24-32%	33-37%	>37%	13-15%	16-23%	24-29%	>29%
over 60	21-24%	25-35%	36-40%	>40%	15-17%	18-26%	27-31%	>31%

Age (years)	Women			Men		
	under-weight/ thin	normal	over-weight/ fat	under-weight/ thin	normal	over-weight/ fat
Athlete 16-50	< 13%	14-26%	>26%	<3%	4-17%	>17%

Guideline for the body water, muscle and bone shares as percentages of the body mass and the calorie requirements to cover the daily requirement:

#### Example for the calculation of the bone percentage:


A man weighs 70 kg and is 35 years old.  
The determined bone mass is  
12.6 kg ( $12.6 \times 100 / 70 = 18\%$ ).

Age (years)	Women				Men			
	Normal values				Normal values			
	Water	Muscles	Bones	Calories (kcal)	Water	Muscles	Bones	Calories (kcal)
6-11	54-65%	24-34%	15-19%	1500-2000	56-67%	25-36%	15-19%	1600-2200
12-15	54-64%	28-38%	14-18%	2000-2400	58-69%	32-44%	15-19%	2200-2700
16-39	50-60%	32-38%	15-19%	2000	55-65%	38-44%	16-20%	2300
40-60	45-55%	30-36%	15-19%	1800	50-60%	36-42%	16-20%	2100
over 60	40-55%	29-35%	14-18%	1600	45-60%	34-40%	15-19%	1800
Athlete 16-50	52-64%	34-41%	15-19%	2600	58-70%	41-48%	16-20%	3000

**Attention:** Owing to great variations in the level of body fat it is difficult to obtain a reference value. On average for adult men the best time for the measurement is between 6 and 8 o'clock p.m. because the level of body fat is very constant during this period!

#### Setting the time and date

After installing the batteries into the display, the time, date, CLOCK and the current ambient temperature will

(fig. 4)

be appear in the display (Fig. 4). Set the time as described below:

Push the **SET** button to set the time. The year will first start to flash on the display. Push either the **P/▲** or **ATH/▼** button. If you hold down the **P/▲** or **ATH/▼** button a quick run will start. Pushing the **SET** button will switch the display to the next setting. Use the

**P/▲** and **ATH/▼** buttons to set the correct month, day, hours and minutes. Confirm your last setting with the **SET** button and the temperature in °C or °F will appear. Pushing the **UNIT** button will change the display from °C to °F and vice versa.

The data must be set again if the batteries are removed from the display or if they must be replaced.

#### Putting into service

##### Simple Weighing Function

Position the balance on a hard, flat surface (avoid carpets). Uneven floors lead to imprecise readings.

1. Switch the scale on by applying slight pressure on the glass surface with your foot or toes. A red indicator lamp lights up thereby signalling that the device is still not ready for use. A green indicator lamp indicates that the scale is ready for use. Switch the display on by pushing the **MODE** button. "8888", WT-FAT, HELO and WT (weight) will appear in the display.

- After several seconds „WT 0.0 kg (0.0 lb./0:0 St:)“ will appear on the display (Fig. 5).

(fig. 5)

- If the scale is not switched on you will see “WT ---- kg (---- lb., ---- St:)” on the LCD display.
- Step on the scale with both feet. Distribute your weight evenly and wait for the scale to calculate your weight. Hold the display vertically or horizontally over the infrared lamp on the scale.
- The reading on the display initially flashes 2x, then your exact weight will be shown. When you step off the scale, the device will switch off automatically after 20 seconds.

## Programming your personal data

**Attention:** Carry out the following operations one after the other and within a few seconds or the device will automatically interrupt the programming mode!

- Push the **MODE** button to switch on the display. Wait until “WT 0.0 kg (0.0 lb./0:0 St:)” appears on the display. Use the **UNIT** button to select your preferred unit (kg, lb. or St:).

- Push the **MODE** button again. Then, push the **SET** button to change to the


(fig. 6)

- programming mode. ‘Set’ appears in the LCD display, then ‘P 01’ (storage location) starts to flash, and the preset symbol for ‘male’ and the preset value for the height in cm (in.) are shown (fig. 6).
- Select a storage location using the **P/▲ / ATH/▼** buttons.
- Then push the **SET** button and the display changes to the next setting.
- The symbol for ‘male’ flashes on the display. Use the **P/▲ / ATH/▼** buttons to

select a setting (‘male’ or ‘female’).

- Push the **SET** button again to change to the HEIGHT setting. The indication ‘170’ (67.0) flashes. Use the **P/▲ / ATH/▼** buttons to select a value between 80 (31.5 in) and 240 cm (94.5 in).
- Push the **SET** button again to change to the AGE setting. The indication ‘30’ flashes. Push one of the buttons **ATH/▼** or **P/▲** and hold it down to run up or down through the values in fast mode. Use the **P/▲ / ATH/▼** buttons to select a value between 6 and 100 years. If the age is set to a value under 16 years, a child’s face will be shown in the display next to the male/female symbol.
- Pushing the **SET** button again will cause the storage location, sex, age and entered height to be displayed for a few seconds, then the indication changes to “WT 0.0 kg (0.0 lb./0:0 St:) and --- (HT)”. If the scale is not switched on the reading will be as shown in Fig. 7.

(fig. 7)

- Push the **MODE** button three times to view all settings again. The reading must be as shown in Fig. 6.
- If you want to make any changes push the **MODE** button until the display reads as shown in Fig. 6. By pushing the **P/▲** button you change to the selection mode. Now choose the storage location with the settings you wish change (flashing indication) using the **P/▲** and **ATH/▼** buttons and push the **SET** button. Proceed as described in steps 4 to 6.
- If you wait for more than 15 seconds without changing anything the reading changes back to the time.

Your personal data is now stored in the selected storage location.


There is an overall number of 16 storage locations. The remaining 15 storage locations can be used to store also the personal data of other members of your family or other users of these scales.

# GB

## Determination of the body and bone mass, and the body fat, muscle and body water percentages

First, enter your personal data in one of the storage locations available. To switch the scale apply pressure on it as described above. Push the **MODE** button to switch on the display; the reading as shown in Fig. 5 will appear. Push the **MODE** button again. By pushing the **P/▲** button you change to the selection mode. Use the **P/▲** / **ATH/▼** buttons to select the storage location where your personal data is stored.

- The display will now first show any stored data as the storage location number (flashing), sex, age (AGE) and body height (HT). Then, "WT 0.0 kg (0.0 lb./0:0 St:) and -.- (HT)" will appear on the display. In this example, the person is 10 years old, female, 137cm tall and occupies storage location 6.

(fig. 8)

### Athlete symbol

This group includes persons who are doing sports for at least 10 hours per week and have a heart rate of about 60 beats/min or lower.

The guidelines do not cover professional sportsmen as, e.g., bodybuilders.

- If the display shows WT 0.0 kg (0.0 lb. / 0:0 St:) and -.- (HT) and your age is set to a value between 16 and 50 years, the **ATHLE** symbol will light up after pushing the **ATH/▼** button.
- Now, step on the scales. Grab the LCD display with both hands and make sure to touch the contact surfaces.
- Wait until the measuring process is completed (the weight indication flashes 1 x).
- After several seconds, the display will show your determined weight, body fat percentage (symbol: "BF %") followed by the muscle mass percentage ("BM %").

Besides this, the scale determines the body water percentage („BW %"), the required calories per day (KCAL) and the bone weight in kg / lb. St: (BON / kg (lb., St:)).

- The determined values will then be shown on the LCD display two times in the sequence described above. Then, the display automatically switches back to the time after some seconds.

### Note:

Only the personal data of a person as, e.g., the age, body height, etc., are stored. The determined values as the body weight, calorie requirement, etc. are only shown once.

## Display and Error Messages

### General

If an error message (Err) appears in the display please wait until the scales switch off. Then switch them on again to continue using them. For a complete reset remove the battery for a moment and then reconnect it.

### Messages

- Switching on: All segments light up simultaneously. This is a check that the displays are working correctly.
- „- - - - kg (- - - - lb. / - - - - St:)": The weight to be measured is too low. Correct operation is not possible. The scale is not switched on.
- „Err": The calculated body fat/water percentage ranges beyond the upper or lower measuring limit. The measurement was aborted before the measuring process was completed or it was not possible to perform measuring.
- „0-Ld": The indication '0-Ld' means that the weight exceeds 150 kg. Do not exceed the device's maximum load to prevent damages to the

device and the sensors.

- Place the scales on a flat and solid surface.

- "Lo-1" The batteries in the display are low.
- "Lo-2" The batteries in the scale are low.

## Usage and care tips

- Weigh yourself without any clothes or shoes on before meals and always at the same hour of the day.
- The balance rounds off weight to the closest 100 g.
- If you weigh yourself twice and two different values are given, it means that your weight is comprised between both values.
- Do not immerse the balance in water: the electronic components can be damaged.
- Remove the batteries if the scales will not be used for a longer period of time.
- Handle the balance with care, as it is a precise measuring instrument. Do not drop it or jump on it.
- Absolutely no need for maintenance.
- Do not load the appliance with more than 150 kg (330,7 lb./23,62 St.) to avoid damaging the device.

## Cleaning

- Clean the scale, i.e. the scale surface and the display, using a slightly moistened, but not wet cloth. Make sure that no moisture or fluid matter gets into the devices!
- Do not use abrasive detergents! These can damage the casing of the device!
- Use a commercially sold detergent or soap for cleaning.
- Thoroughly dry the parts after cleaning using a dry, soft cloth.

This device has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations. Subject to technical changes without prior notice!

## Guarantee

The device supplied by our Company is covered by a 24 month guarantee starting on the date of purchase (receipt).

During the life of the guarantee any fault of the device or its accessories ascribable to material or manufacturing defects will be eliminated free of charge by repairing or, at our discretion, by replacing it. The guarantee services do not entail an extension of the life of the guarantee nor do they give rise to any right to a new guarantee!

Proof of the guarantee is provided by the proof of purchase. Without proof of purchase no free replacement or repair will be carried out.

If you wish to make a claim under the guarantee please return the entire machine in the original packaging to your dealer together with the receipt.

Defects to consumables or parts subject to wearing as well as cleaning, maintenance or the replacement of said parts are not covered by the guarantee and hence are to be paid!

The guarantee lapses in case of unauthorized tampering.

## After the expiry of the guarantee

After the expiry of the guarantee repairs can be carried by the competent dealer or repair service against the payment of the ensuing costs.

## Meaning of the “Dustbin” Symbol

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Use the collection points provided for the disposal of electric appliances to dispose of electric appliances that you do not use any more.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of reutilisation of electrical and electronic equipment.  
Information concerning where the equipment can be disposed of can be obtained from your local authority.

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